



TECHNOPATH
CLINICAL DIAGNOSTICS

QUALITY CONTROLS FOR OPTIMAL PATIENT CARE



YOUR QUALITY CONTROL SOFTWARE

A background image showing a man and a woman in white lab coats looking at a laptop screen. The image is overlaid with a blue digital network pattern, including glowing nodes, connecting lines, and floating numbers (1-9).

EXPERT

CONTENTS

IAMQC Introduction	3
Work Flow	4
Technopath Products	5
Expert	6
Connectivity	10
Benefits	11






Designed to complement and support TECHNOPATH's Multichem™ Quality Control (QC) product range, IAMQC™ Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time.

IAMQC™ Software products allow users to automate, centralise, standardise and improve QC processes in a laboratory setting. Our combination of modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

TECHNOPATH's full suite of software products provide clinical laboratories significant cost and time savings, whilst delivering higher confidence in analytical testing methods. Choose from Intranet and/or Internet based statistical quality control and quality assurance software products. IAMQC™ software products are practical, graphical, user-definable and easy to use.





TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH'S SOFTWARE PRODUCTS HELP:

BENCH TECHNOLOGISTS:


- Spend less time on false positive QC flags
- Concentrate on tests, which require their attention
- Spend less time trouble-shooting
- Know how to react when the mean shifts
- Assess the acceptability of new reagent lots and calibrations
- Solve QC problems
- Gain understanding and confidence in the QC process

LAB MANAGERS:

- Choose QC rules to maximise true rejects and minimise false rejects
- Quickly see the tests that require their attention
- Skim graphics to quickly review current or historical data by lab, department, instrument or test
- Monitor performance in groups of laboratories
- Review problem tests and QA activities in local and remote labs

LAB OR HOSPITAL ADMINISTRATORS:

- Save money
- Improve quality
- Improve service
- Review Administrative Summary Reports to ensure quality performance



TECHNOPATH
CLINICAL DIAGNOSTICS
QUALITY CONTROLS FOR OPTIMAL PATIENT CARE


EXPERT

IAMQC™ Expert is an interactive system that helps front-line laboratory staff select QC rules, reduce unnecessary repeats and make meaningful QC decisions. IAMQC™ Expert creates graphical color-coded representations of accuracy and precision, illustrating performance relative to target values and error limits. Interactive modules such as 'Reagent Verification', 'Calibration Check' and 'Mean Shift Analysis' assist with problem solving and decision-making.

IAMQC™ Expert allows the end-user to monitor method performance relative to clinical requirements and focus on the tests that require their attention. The system can also be integrated with IAMQC™ Daily to automatically capture summary data for analysis. Upon entering the IAMQC™ Expert database, the software automatically analyses the data and suggests QC rules to maximise true rejects and minimise false rejects.

KEY FEATURES

- Works with multiple Sites, Departments, Instruments, Tests and Levels.
- Recommends QC rules for IAMQC™ Daily software.
- Powerful QC troubleshooting tools.
- Integrates with various Laboratory Information Systems and/or instruments.
- Works on a single PC, LAN, WAN, and over the Internet.
- Runs on a powerful Database Management System.
- Powerful reporting and charting capabilities




EXPERT ANALYSIS

Times and technologies are changing rapidly. Instruments and methodologies are more accurate, precise and stable than they were a decade ago. Most laboratories have adopted these new technical advances, but few have modified their QC processes to match. Many laboratories are still using a 1-2s rule as recommended by Levey and Jennings in 1951.

In 1981 Westgard recommended using a multi-rule algorithm to avoid the false positive flags inherent in using only a 1-2s rule. Technology has come a long way since 1981, and for the past ten years industry leaders such as James Westgard, Per Hytøft Petersen and Callum Fraser have advised the use of a variety of QC rules to match the analytical capability and stability of each test. With the dramatically improved precision of today's methods, we now see shifts of several SD for the same control on the same test from time to time within a single laboratory. We have designed a QC system that will alert users to significant changes and not generate QC flags when the system is operating safely within acceptable limits.

The system compares method performance to defined quality requirements (rather than to last month's data) and recommends QC strategies that will warn users when QC data points exceed acceptable performance - with a minimal number of false flags. In the design of our QC system we "balance" the quality control system to meet the changing performance and stability of the analytical system.

Our analytical processes also vary in their susceptibility to the occurrence of significant errors. Error rates vary from low to moderate to high, depending on the frequency of significant errors that occur in a specific test system. Some methods seldom encounter significant problems. Others are susceptible to relatively frequent sources of significant error. Methods with high error rates frequently see significant shifts in the mean or have on-going precision problems; these methods may also be susceptible to frequent instrument breakdowns or problems.



Data Review - Traffic Light System

We monitor method performance (accuracy and precision) relative to a quality requirement by calculating critical systematic error (\wedge SEc). Critical systematic error is an extremely powerful and useful statistic. \wedge SEc indicates in one number how method accuracy and precision compare to the target and TEa limit set for each control. \wedge SEc indicates the number of standard deviations the mean can shift before the results will exceed error limits. Therefore changes in either the mean or SD will be reflected in a change in critical systematic error.

Our QC system has an improved error detection process when methods are close to the error limit (have a low \wedge SEc) or when methods have poor stability and are prone to errors. In this case we will select a higher number of controls or run our existing controls more frequently and we will select QC rules that are more powerful when assessing small changes in method performance. When the \wedge SEc is high (method performance is well within quality requirements) and the analytical system is stable, we will run fewer controls or run controls less frequently

and we will select QC rules to minimise false QC flags.

Critical systematic error is a valuable indicator of the size of the shift in the mean that we must detect. We can easily visualize how it would be appropriate to use "tighter" QC rules when the method is closer to the error limit and to use "looser" QC rules when a method can shift many standard deviations before exceeding the quality requirement. When we run more levels of controls or run the same controls more frequently, we increase our probability of detecting errors. Unfortunately, we also increase the probability of false rejection. Remember, if we are using a 1-2S rule, we will see 5% of our "good" data falling between 2 and 3 SD. Therefore, the more controls we run, the higher the probability in any given run that one of them will fall outside two standard deviations. When it is necessary to design a QC process to detect a very small change in the analytical system, one of the strategies we can use is to increase the number and frequency of the controls.

Microsoft Excel - Problem Solver Single1

File	Edit	View	Insert	Format	Tools	Data	Window	Help																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																					
------	------	------	--------	--------	-------	------	--------	------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

IAMQC™ Expert is a user-friendly interactive expert system that helps frontline laboratory staff select QC rules, reduce unnecessary repeats and make meaningful QC decisions. The software creates graphical colour-coded representations of accuracy and precision illustrating performance relative to target values and error limits for several departments or groups of affiliated laboratories. Interactive modules assist with problem solving and decision-making.

IAMQC™ Expert Helps Bench Technologists:

- Spend less time on false positive QC flags
- Concentrate on tests which require their attention
- Spend less time trouble-shooting
- Know how to react when the mean shifts
- Assess the acceptability of new reagent lots and calibrations
- Solve QC problems
- Gain understanding & confidence in the QC process

IAMQC™ Expert Helps Lab Managers:

- Monitor method performance relative to CLIA or clinical requirements
- Choose QC rules to maximise true rejects and minimise false rejects
- Quickly see the tests that require their attention
- Skim graphics to quickly review current or historical data by lab, department, instrument or test
- Monitor performance in groups of laboratories
- Review problem tests and QA activities in local and remote labs

IAMQC™ Expert Helps Lab Or Hospital Administrators:

- Save money
- Improve quality
- Improve service
- Review Administrative Summary Reports to ensure quality performance

Microsoft Excel - DetailQualContDataLab1

FileEditViewInsertFormatToolsDataWindowHelp

Type a question for help

110%

Times New Roman

14

B

I

U

</

Expert Reports in Microsoft Excel

CONNECTIVITY

Providing the flexibility for an internet OR intranet connection facilitates the customer in choosing their preference with regards to connectivity. The system can run locally behind a firewall, or alternatively, over the web should access be required outside of the local network.

Built-in real-time and semi-real time interface solutions are available to capture data from all types of instruments,

middleware systems and LIS across all departments. To date, we have developed over 200 types of drivers for data capturing purposes. IAMQC™ Daily can also capture QC results from diagnostic instruments that are not interfaced to the DMS and from manual result entry programs on a daily basis. IAMQC™ Drivers include, but are not limited to, the following list:

Abbott Architect	Dawning UC	RCM Beziers
Abbott AXSYM	EPIC LIS SDF 1.0	Roche COBAS
Abbott CELL DYN	ERMA	Roche ELECSYS
ABX	Fletcher Flora	Roche HITACHI 747
AVL OMNI 5	HMS LIS	Roche Integra 4.1.ARC
Bayer Atlas	Hutt No 2 QC	Roche Modular
Bayer CLINITEK	IL ELECTRA 1000	Roche PSM
Beckman AU	IL Synthesis	RT Communicator
Beckman DxC	Immucor Galileo	Schuytlab LIS
Beckman DXI ACCESS 2	LabDaq LIS	Siemens ADVIA Centaur
Beckman Remisol	McKesson LIS	Siemens Centralink
Beckman SYNCHRON CX9	Meditech LIS	Siemens Dimension 2.0
Cerner LIS	MGC 240 Qualitative 1.0	Siemens Immulite (DPC)
Cerner Millennium LIS	Microplate	Siemens Novius LIS
Citation LIS	Microscan	SOFT LIS
CLARIS LIS	MOLIS LIS	StaRRsed
Coag-A-Mate MTX II	Mysis LIS	Stat Profile M
ConcurTrak SDF	NEMO Middleware	Sunquest 1
CPSI LIS	NOVA ELECTROLYTE 5	SYMEX
Dade Behring BN II Nephelometer	NUCLEUS	TOSOH A1c 2.2
Dade Dimension	Omnilab AMS	UF100
Dairyland LIS	Orchard LIS	UriScan-S300
Data Innovations IM	Ortho Vitros Eci	Viper
Datalink	Quadrated LIS	RCM Beziers

BENEFITS

- ① **Centralised review** of all QC data from all laboratories/ instruments. Central administrator access to review QC performance at multiple facilities – no need to visit each laboratory site
- ② **Closer monitoring** of QC from remote locations without additional costs
- ③ **Built-in real-time and semi-real time interface solutions.** Integrate with various Laboratory Information Systems and/or instruments
- ④ **Capture QC results** from diagnostic instruments and manual result entry programs on a daily basis
- ⑤ Compare each result to **Assigned Mean & SD**
- ⑥ **Assess QC results** against a set mean and standard deviation using QC rules (Westgard and/or User-defined): single or multiple QC rules
- ⑦ Auto-approval protocol
- ⑧ Troubleshoot problematic daily QC
- ⑨ **Manage multiple Sites, Departments, Instruments, Tests and QC Levels on one central database**
- ⑩ **Manage both quantitative and qualitative results**
- ⑪ **Manage different departments** (Chemistry, Haematology, Microbiology, etc.) on one software system
- ⑫ **Focused troubleshooting** for failed QC results
- ⑬ **Technologist and Supervisor Reviews/Sign-off**
- ⑭ **"Reverse Levels"** automatic function
- ⑮ **QC management at different levels:** administrative and bench technologists
- ⑯ **Document all activities** regarding daily QC
- ⑰ **Tracking of proficiency testing performance and problem resolution**
- ⑱ **Documentation** of new reagent/calibrator/QC lot numbers and studies
- ⑲ **Monthly reporting on-line** for management
- ⑳ **Document activities and administrative comments** for summarized QC data
- ㉑ **Transfer detail and summary data** between a Laboratory and a single QC database in real time over the Internet
- ㉒ **Internal/External peer QC review capability.** Collect, analyse and compare individual laboratory data immediately with a world-wide peer group at the touch of a button over the Internet
- ㉓ **Works on a single PC, LAN, WAN, and over the Internet**
- ㉔ **Runs on a powerful Database Management System** to support large volumes of data in real time
- ㉕ **Multiple Assign and Advanced Setup/Copy functions** for fast and easy setup and ongoing maintenance
- ㉖ **Multiple ways to enter data manually** (by Level, Test or Instrument; one at a time or many at a time)
- ㉗ **Monthly Supervisor Review**
- ㉘ **Powerful reporting and charting capabilities.** Includes User-definable reports
- ㉙ **Multi-user environment**
- ㉚ **Audit trail/Admin module** for setting up users with different security profiles





TECHNOPATH
CLINICAL DIAGNOSTICS

Fort Henry Business Park, Ballina, Co. Tipperary, Ireland.
info@technopathCD.com www.technopathclinicaldiagnostics.com