# Laboratory Evaluation of Technopath Controls on the **ARCHITECT Chemistry Analyzers**

D. Brault<sup>1</sup>, L. Lennartz<sup>2</sup>, M. Orth<sup>3</sup>, J. Shih<sup>4</sup>, and A. Croce<sup>5</sup>

<sup>1</sup>Hôpital Tenon, Service de Biochemie, Paris, France; <sup>2</sup>Abbott Laboratories, Wiesbaden, Germany; <sup>3</sup>Instituts for Laboratoriumsmedizin, Vinzenz von Paul Klinikem, Marienhospital, Stuttgart, Germany; <sup>4</sup>Abbott Laboratories, Abbott Park, IL, US; <sup>5</sup>Ospedale Civile Sondrio, Sondrio, Italy

2014 IFCC WorldLab Meeting, Istanbul • 22 - 26 June

#### **Background**

Efficient analytical and quality control with the BioRad Liquichek and BioRad ABBOTT solution is important for testing laboratories to ensure released results meet the required quality in regard to accuracy and precision. Technopath Manufacturing Ltd. recently multi-constituent introduced controls (MCC) for use with the Abbott ARCHITECT instruments. The performance of the urine and serum chemistry controls were evaluated in comparison to the laboratory's current QC.

#### **Materials and Methods**

Technopath Multichem S Plus controls are prepared from human serum to which purified biochemical material (extracts of human and animal origin), chemicals, drugs, preservatives and stabilizers are added. These controls are provided in liquid form and stored frozen (-20 to -80°C) until use. Once thawed, the material is stored at 2–8°C and stable for 10 days unless otherwise stated in the lot specific data sheets. Sixty analytes are included in this control at three different levels.

Multichem U controls are prepared from human urine to which purified biochemical material (extracts of human and animal orign), chemicals, drugs, preservatives and stabiliers have been added. These controls are provided in liquid form and stored at 2–8°C. Once the material is opened, it should be stored tightly capped at 2-8°C and is stable for 30 days unless otherwise stated in the lot specific data sheets. Thirteen analystes are included in this control at two different levels.

Three European sites (Paris, France; Multichem S Plus Stuttgart, Germany; Sondrio, Italy) evaluated the Technopath Multichem S Plus and Multichem U controls for a minimum of thirty days in parallel with the lab's routine QC controls including single constituent controls (SCC) and MCC for assays on their routine menus. Testing was performed with both the ARCHITECT c8000 and c16000 instruments.

Data presented here are from the following serum analytes: AALT\*, AAST\*, total bilirubin, Cl, total cholesterol, glucose, K, total protein, Na, triglycerides, and urea; and the following urine analytes: Cl, creatinine, glucose, K, Na, and urea. All data were collected via AbbottLink for automated data retrieval. Means, standard deviation and range were calculated for all controls. Assay reagent lots and calibrator lots varied across the sites and within the sites.

Other controls evaluated included BioRad Liquichek Unassayed Chemistry controls, Liquichek Lipids, Liquichek Ethanol/ Ammonia, Liquid Assay Multiqual, and Liquichek Urine Chemistry.

Sigma Metrics were calculated using peer means and CVs and TEas from the literature according to the following equation: Sigma Multichem Urine Metric = (TEa - Bias)/CV.

#### Results

Results from 12 frequently performed clinical chemistry assays were analyzed. The %CV for the 12 assays with the Multichem S Plus control ranged from 0.42 (chloride) to 4.71% (total bilirubin).

The %CV for the 6 assays with the Multichem U control ranged from 0.50 (chloride and potassium) to 5.24% (Urea).

For both controls, the majority of the CVs were less than 2%. The results from the Technopath controls compared favorably

### Results (cont'd)

Multiqual serum controls and the BioRad Liquichek Urine control in terms of %CV.

Little variation was seen instrument to instrument or even site to site. When all the data was consolidated by analyte across multiple reagent lots and instruments, the overall %CV ranged from 0.765 (chloride) to 4.11% (AALT) for S Plus and 1.01 (sodium) to 3.171% (glucose) for Urine.

#### **Summary of Sites and Testing -Multichem S Plus**

	France	Italy	Germany
Instruments	2 - c8000	2 - c8000 1 - c16000	2 - c8000
AALT	Χ		X
AAST	X		X
Bilirubin, Total	X	X	X
Chloride	Χ	X	X
Cholesterol, Total	Χ	Χ	X
Creatinine, Enzymatic	X		
Creatinine, Picrate		Χ	X
Glucose	X	X	X
Potassium	Χ	Χ	X
Protein, Total	X	X	X
Sodium	Χ	X	X
Triglycerides	Χ	X	X
Urea	Χ	X	Χ

#### Summary of Sites and Testing – **Multichem U**

	France	Italy	Germany
Instruments	2 - c8000	2 - c8000 1 - c16000	2 - c8000
Chloride	Χ	Χ	Χ
Creatinine, Enzymatic	Χ		
Creatinine, Picrate		X	X
Glucose	Χ	X	Χ
Potassium	Χ	Χ	Χ
Sodium	Χ	X	Χ
Urea	Χ	Χ	Χ

Analyte	Level	Unit	N	Expected Mean	Mean	SD	cv	Expected Range	Observed Range
AALT	1	U/L	153	26.9	27.1	1.12	4.11	21.5-32.3	24.0-30.0
	2		156	114	114.9	3.18	2.77	91.2 - 137	104.7 - 119.0
	3		153	250	250.2	5.35	2.14	200-300	235.1-259.0
AAST	1	U/L	149	45	44.7	1.8	4.04	36.0-54.0	40.0-50.0
	2		150	132	131.9	3.35	2.54	106-158	122.0-142.0
	3		149	261	264.7	4.52	1.71	209-313	255.0-281.0
Bilirubin	1	µmol/L	267	17.2	16.98	0.684	4.03	13.8 – 20.6	14.54 – 19.15
	2		266	45.6	44.8	1.284	2.87	36.5 – 54.7	39.84 – 48.22
	3		259	95.7	91.89	3.611	3.93	76.6 – 115	80.37 – 98.15
Chloride	1	mmol/L	260	85.1	83.7	0.904	1.08	76.6 – 93.6	81.50 – 88.20
	2		260	95.6	93.59	0.806	0.86	86.0 – 105	90.70 – 96.50
	3		255	110	108.22	0.828	0.76	99.0 – 121	105.70 – 111.60
Cholesterol	1	mmol/L	248	2.66	2.7	0.042	1.54	2.13 – 3.19	2.50 - 2.83
	2		244	4.06	4.09	0.087	2.12	3.25 - 4.87	3.61 – 4.37
	3		238	6.43	6.49	0.109	1.68	5.14 – 7.72	6.14 - 6.73
Creatinine, E	1	µmol/L	85	69.4	69.2	2.11	3.05	55.5 - 83.3	62.7 – 72.5
	2		90	175	176.2	3.18	1.8	140 – 210	167.6 – 183.8
	3		89	519	525.2	6.35	1.21	415 – 623	512.6 - 551.7
Creatinine, P	1	mg/dL	148	0.697	0.7	0.023	3.32	0.558 <b>–</b> 0.837	0.64 - 0.79
	2		147	1.97	1.98	0.062	3.12	1.57 - 2.36	1.69 - 2.26
	3		145	5.02	6.13	0.11	1.79	4.82 - 7.23	5.92 - 6.63
Glucose	1	mmol/L	259	2.82	2.84	0.049	1.72	2.26 - 3.38	2.70 - 3.02
	2		258	7.37	7.4	0.134	1.81	5.90 - 8.84	6.93 - 7.89
	3		253	16.2	16.21	0.236	1.46	13.0 - 19.4	15.46 - 17.05
Potassium	1	mmol/L	262	2.42	2.45	0.051	2.07	2.18 - 2.66	2.30 - 2.59
	2		260	3.77	3.8	0.058	1.52	3.39 – 4.15	3.69 - 3.94
	3		257	6.76	6.75	0.072	1.07	6.08 - 7.44	6.45 - 7.06
Protein, T	1	g/L	190	46.1	45.66	0.961	2.11	36.9 – 55.3	41.70 - 48.30
	2		190	65.2	64.9	1.329	2.05	52.2 – 78.2	60.10 - 68.70
	3		189	87.5	87.83	1.067	1.21	70.0 – 105	81.40 - 92.60
Sodium	1	mmol/L	259	121	120.68	1.068	0.88	109 – 133	117.60 – 125.00
	2		259	141	139.88	1.199	0.86	127 – 155	134.80 – 144.50
	3		256	161	160.39	1.264	0.79	145 – 177	156.90 – 165.00
Triglycerides	1	mmol/L	194	0.625	0.63	0.019	2.98	0.500 - 0.750	0.56 - 0.68
	2		190	1.5	1.47	0.037	2.53	1.20 - 1.80	1.36 - 1.55
	3		188	2.52	2.48	0.046	1.84	2.02 - 3.02	2.30 - 2.59
Urea	1	mmol/L	78	3.05	3	0.09	3.05	2.44 - 3.66	2.8 - 3.2
	2		80	14	13.7	0.29	2.11	11.2 - 16.8	12.7 - 14.8
	3		78	22.8	22.3	0.49	2.19	18.2 - 27.4	20.7 - 23.7

Analyte	Level	Unit	N	Expected Mean	Mean	SD	CV	Expected Range	Observed Range	
Chloride	1	mmol/L	263	87.2	87.9	1.59	1.8	78.5 – 95.9	78.0 – 96.0	
	2		256	193	189.4	2.22	1.17	174 – 212	181.0 – 201.0	
Creatinine, E	1	mmol/L	89	6.42	6.05	0.11	1.81	5.14 – 7.70	5.73 - 6.37	
	2		88	12.9	11.99	0.205	1.71	10.3 – 15.5	11.51 – 12.72	
Creatinine, P	1	mg/dL	72	66.1	66.06	1.598	2.42	52.8 – 79.3	63.31 – 72.02	
	2		71	135	131.63	2.302	1.75	108 – 162	125.60 – 136.63	
Glucose	1	mmol/L	262	1.99	2.03	0.064	3.17	1.59 – 2.39	1.71 – 2.23	
	2		264	20	20.41	0.543	2.66	16.0 – 24.0	17.38 – 21.69	
Potassium	1	mmol/L	259	16.2	16.2	0.31	1.93	14.6 – 17.8	13.8 – 17.7	
	2		259	58.8	58.9	0.82	1.4	52.9 – 64.7	54.8 - 61.9	
Sodium	1	mmol/L	259	81.8	81.5	1.7	2.09	73.6 – 90.0	76.0 – 89.0	
	2		259	158	157.9	1.6	1.01	142 – 174	151.0 – 168.0	
Urea	1	mmol/L	79	153	151.4	2.97	1.96	122 – 184	145.3 – 162.4	
	2		79	358	347.4	7.59	2.18	286 – 430	329.4 – 371 6	

#### ARCHITECT MCC Study - Assay: Activated Alanine **Aminotransferase (AALT) (8L92)**





Figure 1: Box and Whisker plots of activated quality control testing for routine controls - BioRad assayed and unassayed Multiconstituent Controls (MCC) - and Technopath S Plus

#### **ABBOTT ARCHITECT Technopath MCC** Study - Assay: Chloride (S) (2P32)



Box and Whisker plots of Chloride quality control testing for routine controls - BioRad assayed and unassayed Multiconstituent Controls (MCC) - and Technopath S Plus

#### **ABBOTT ARCHITECT Technopath MCC** Study - Assay: Potassium (U) (2P32)




Figure 3: Box and Whisker plots of Urine Potassium quality control testing for routine controls - BioRad Liquichek Urine Chemistry and Technopath U

#### **ABBOTT ARCHITECT Technopath MCC** Study - Assay: Sodium (U) (2P32)




Figure 4: Box and Whisker plots of Urine Sodium quality control testing for routine controls - BioRad Liquichek urine chemistry controls - and Technopath Multichem U control

#### Technopath AALT Method **Decision Chart**

# **Decision Chart**


**Chloride Method** 

Figure 5. Six Sigma Figure 6. Six Sigma Method Performance Method performance comparison of BioRad comparison of BioRad (o) and Technopath (o) and Technopath (o) ( controls for AALT. controls for chloride.

#### **Sigma Metric Summary: Multichem S Plus**

	Olgille	i ivicui	c – iec	Jaui	Olgina Wetric - Dioriau					
Analyte	6	5	4	3	<3	6	5	4	3	<3
AALT	15	3	2		1	6	3		1	4
AAST	19		2			12	1	1		
Bilirubin, Total	14	4	3			9	4	1		
Chloride	12	5	1	3		5	3	2	4	
Cholesterol	7	3	5	3		9	2	2		
Creatinine E	5			1		3	1			
Creatinine P	10		1	1		7	1			
Glucose	13	5	1	2		9	2	1	2	
Potassium	18	2	1			10	1	1	2	
Sodium		6	3	5	7			3	7	4
Total Protein	10		3	1	1	8		1		
Triglycerides	15					10				
Urea	1	2	4	5	9		1	3	3	7
Total	139	30	26	21	18	88	19	15	19	15
%	59.4	12.8	11.1	9	7.7	56.4	12.2	9.6	12.2	9.6

#### **Urine Potassium Method Decision** Chart




( controls. For urine ( controls. For urine Potassium.

Chart Urine Sodium Control Comparison, RCPA10% Goa

**Method Decision** 

**Urine Sodium** 



## Sigma Metric Summary: Multichem U

Sodium.

	Sigma Metric - Technopath						Metric – Technopath   Sigma Metric – BioRad				
Analyte	6	5	4	3	<3	6	5	4	3	<3	
Chloride	11	1	1	1	0	13	1		0	0	
Creatinine E	2		2			2	1	1			
Creatinine P	3	1	2		2	3	0	2	2	1	
Glucose	8	1	2	1	2	10	4	0	0		
Potassium	11	1	0	2		10	4	0			
Sodium	10		3	0	1	11	0	1	1	1	
Urea	12	1	0	1		10	2	1	1		
Total	57	5	10	5	5	59	12	5	4	2	
%	69.5	6.1	12.2	6.1	6.1	72	14.6	6.1	4.9	2.4	

#### **Conclusions**

The Technopath Multichem S Plus and U controls showed comparable performance to the routine QC solution at the three sites.

The new MCCs allow the optimization of day-to-day analytical QC testing by having routine metabolites, specific proteins and TDMs consolidated into one control.

ABBOTT ARCHITECT c8000 and c16000 are trademarks of Abbott Laboratories.

<sup>\*</sup> These assays are not available in the U.S.