



TECHNOPATH

CLINICAL DIAGNOSTICS

Multichem[®] CSF

Providing Third-Party Cerebral Spinal Fluid QC
in a Liquid Stable Format



Features

- ◆ Human based proprietary formulation
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 36 month closed vial stability at -20°C to -80°C .
- ◆ 30 day open vial stability at 2°C to 8°C .

Product	Description	Configuration	Part Code
Multichem CSF	Unassayed Bi-Level	2 × 6 × 2mL	CF100X
	Assayed Bi-Level	2 × 6 × 2mL	CF100A

Multichem[®] CSF Analyte List

Analytes

Glucose
Lactate
IgG
Protein



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathcd.com



TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS

IRELAND

Technopath Life Sciences Park,
Fort Henry, Ballina, Co.Tipperary, V94 FF1P

USA

99 Lafayette Drive,
Suite 179 Syosset, NY 11791

For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathcd.com