



TECHNOPATH

CLINICAL DIAGNOSTICS

Multichem[®] A1c

Providing Third-Party Diabetes Haemoglobin A1c QC in a Liquid Stable Format



Features

- ◆ Human based proprietary formulation.
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 24 month closed vial stability at -20°C to -80°C.
- ◆ 30 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem A1c	Assayed Single Level (Level 1)	12 × 1mL	HB001A
	Assayed Single Level (Level 2)	12 × 1mL	HB002A
	Assayed Bi-level	2 x 6 x 1mL	HB000A

Multichem[®] A1c Analyte List

Analytes

A1c

% HbA1c



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathcd.com



TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS

IRELAND

Technopath Life Sciences Park,
Fort Henry, Ballina, Co. Tipperary, V94 FF1P

USA

99 Lafayette Drive,
Suite 179 Syosset, NY 11791

For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathcd.com