

TECHNOPATH CLINICAL DIAGNOSTICS

Multichem[®] **P**

Providing Third-Party Supplementary Immunoprotein QC in a Liquid Stable Format

TECHNOPATH Multichem[®] P





3ml TECHNOPATH, Ballina, Co. Tipperary

Features

- Human based proprietary formulation.
- Mimics performance of patient samples.
- Targeted at key clinical decision points.
- Liquid for ease of use.

Specifications

- ♦ 36 month closed vial stability at -20°C to -80°C.
- ♦ 14 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem P	Unassayed Single Level	12 × 3mL	SP40PX
	Assayed Single Level	12 × 3mL	SP41PA

THE QUALITY CONTROL COMPANY

Multichem[®] P Immunoprotein Analyte List

Immunoproteins

Alpha-1 Acidglycoprotein Alpha-1 Antitrypsin Alpha-2-Macroglobulin* Antistreptolysin O (ASO)* ADNase B (Anti-Streptococcal DNase B)* Antithrombin III* Apolipoprotein A1 (APO A1) Apolipoprotein B (APO B) Beta-2 Microglobulin C1 Inhibitor* CH50 (Total hemolytic Complement)* Cystatin C*

Chemistry Analytes

Albumin* Angiotensin Converting Enzyme* Total Protein* Complement C3 Complement C4 Ceruloplasmin C-Reactive Protein Ferritin* Haptoglobin Hemopexin* Immunoglobulin A Immunoglobulin G Immunoglobulin M

lgE*

IgG1, Subclass* IgG2, Subclass* IgG3, Subclass* IgG4, Subclass* Kappa Light Chain* Lamda Light Chain* Lipoprotein (a)* Prealbumin Properdin Factor B* Retinol Binding Protein* Rheumatoid Factor Transferrin sTIR (Soluble Transferrin Receptor)*

*Please refer to lot specific package inserts for stability and performance claims.



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in realtime. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathcd.com



TECHNOPATH CLINICAL DIAGNOSTICS Technopath Life Sciences Park, Fort Henry, Ballina, Co.Tipperary, V94 FF1P, Ireland. For Multichem QC Support: qcsupport@technopathcd.com Website: www.technopathcd.com