



TECHNOPATH
CLINICAL DIAGNOSTICS
USA

Multichem[®]

IA Speciality

Providing Third-Party Specialty Peptide
Hormone QC in a Liquid Stable Format



Features

- ◆ Human based proprietary formulation
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 36 month closed vial stability at -20°C to -80°C .
- ◆ 7 day open vial stability at 2°C to 8°C .

Product	Description	Configuration	Part Code
Multichem IA Speciality	Unassayed Single Level (Level 1) Multichem IA Speciality Control	12 x 2mL	BP301X
	Unassayed Single Level (Level 2) Multichem IA Speciality Control	12 x 2mL	BP302X
	Unassayed Single Level (Level 3) Multichem IA Speciality Control	12 x 2mL	BP303X
	Assayed Single Level (Level 1) Multichem IA Speciality Control	12 x 2mL	BP301A
	Assayed Single Level (Level 2) Multichem IA Speciality Control	12 x 2mL	BP302A
	Assayed Single Level (Level 3) Multichem IA Speciality Control	12 x 2mL	BP303A
Unassayed Tri-Level Multichem IA Speciality Control	Unassayed Tri-Level Multichem IA Speciality Control	3 x 4 x 2mL	BP300X
	Assayed Tri-Level Multichem IA Speciality Control	3 x 4 x 2mL	BP300A
Dropper Tips	Optional Dropper Tip	1 x 100pcs	987000

THE QUALITY CONTROL COMPANY

Multichem® IA Speciality

Analyte List

Analytes

ACTH
BNP
Calcitonin
Procalcitonin
PTH Intact



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyze their QC results in real-time. IAMQC Software products allow users to automate, centralize and standardize QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathusa.com



TECHNOPATH
CLINICAL DIAGNOSTICS
USA



TECHNOPATH CLINICAL DIAGNOSTICS USA
99 Lafayette Drive, Suite 179 Syosset, NY 11791
For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathusa.com

JV_MM_IASpeciality_Digi_R3