



TECHNOPATH
CLINICAL DIAGNOSTICS
USA

Multichem® ID

INFECTIOUS DISEASE QC

POWERED BY



Multichem® ID-B

Providing Instrument Specific Third-Party
Infectious Disease Quality Controls



Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control tubes
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

Specifications

- ◆ 24 month closed tube stability at 2°C to 8°C.
- ◆ 90 day open tube stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-B	Positive Control	4 x 5mL	SR101B
Multichem ID-B Mini Kit	Positive Control Mini Kit	1 x 5mL	SR101MB

For further information see Instructions For Use.

INFORMATION FOR USA ONLY:
This product is not intended for use
in donor screening testing.

THE QUALITY CONTROL COMPANY

Multichem® ID-B

Analyte List

Analytes

Anti-HIV 1 IgG
 Anti-HCV IgG
 Anti-HBc IgG
 HBsAg
 Anti-HTLV IgG
 Anti-Treponema Pallidum IgG

Analyser Systems

MANUFACTURER	SYSTEM	ASSAY NAME
Abbott	ALINITY i	Anti-HBc II
		Anti-HCV
		HBsAg Qualitative II
		HIV Ag/Ab Combo
		rHTLV I/II
		Syphilis TP
Abbott	ALINITY s	Anti-HBc
		Anti-HCV
		HBsAg
		HIV Ag/Ab Combo
		HTLV I/II
		Syphilis
Abbott	ARCHITECT	Anti-HBc II CMIA
		Anti-HCV CMIA
		HBsAg Qualitative II CMIA
		HIV Ag/Ab Combo CMIA
		rHTLV-I/II CMIA
		Syphilis TP CMIA
DiaSorin	LIAISON	Anti-HBc CLIA
		Treponema Screen CLIA
DiaSorin	LIAISON XL MUREX	HBsAg Quant CLIA
		HCV Ab CLIA
		HIV Ab/Ag CLIA
		recHTLV-I/II CLIA



TECHNOPATH
 CLINICAL DIAGNOSTICS
 USA

TECHNOPATH CLINICAL DIAGNOSTICS USA
 99 Lafayette Drive, Suite 179 Syosset, NY 11791
 Website: www.technopathcd.com