



TECHNOPATH
 CLINICAL DIAGNOSTICS
 USA

Multichem® ID

INFECTIOUS DISEASE QC

POWERED BY



Multichem® ID-SeroNeg

Providing Instrument Specific Third-Party
 Infectious Disease Quality Controls



Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control tubes
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

Specifications

- ◆ 36 month closed tube stability at 2°C to 8°C.
- ◆ 90 day open tube stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-SeroNeg	Negative Control	4 x 5mL	SR100N
Multichem ID-SeroNeg Mini Kit	Negative Control Mini Kit	1 x 5mL	SR100MN

For further information see Instructions For Use.

INFORMATION FOR USA ONLY:
 This product is not intended for use
 in donor screening testing.

Multichem[®] ID-SeroNeg

Analyte List

Analytes

Anti-HIV

Anti-HCV

Anti-HBc

Anti-HTLV

Anti-Treponema

HBsAg

HIV p24



TECHNOPATH
CLINICAL DIAGNOSTICS
USA

TECHNOPATH CLINICAL DIAGNOSTICS USA
99 Lafayette Drive, Suite 179 Syosset, NY 11791
Website: www.technopathcd.com

JV_MM_IDSeroNeg_Digi_R1