



**TECHNOPATH**  
CLINICAL DIAGNOSTICS  
USA

# Multichem<sup>®</sup> ID-COVID19 G

Unassayed Third-Party Quality Control for  
Antibodies to SARS-CoV-2 (including IgG)

## Features

- Ready to use liquid stable controls
- Technical and scientific support

## Specifications

- 24 month closed vial stability for CVG200P\* at 2°C to 8°C.
- 21 month closed vial stability for CVG300P\* at 2°C to 8°C.
- 21 day open vial stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-COVID19 G	Unassayed Positive Control	4 x 4mL	CVG200P
			CVG300P
Multichem ID-COVID19 G	Unassayed Positive Control	4 x 2mL	CVG100P

## Analyte

Positive for Antibodies to SARS-CoV-2 (including IgG)

For further information see Instructions For Use.

Note: Additional information on positive QC for SARS-CoV-2 IgM is available by contacting [multichemcovid@technopathcd.com](mailto:multichemcovid@technopathcd.com)

\* Lot Specific

*THE* QUALITY CONTROL COMPANY

**SEROLOGY  
CONTROLS**

**Multichem<sup>®</sup> ID**

INFECTIOUS DISEASE QC



**TECHNOPATH CLINICAL DIAGNOSTICS**

[www.technopathcd.com](http://www.technopathcd.com)

[multichemcovid@technopathcd.com](mailto:multichemcovid@technopathcd.com)

MM\_IDCOVID19\_G\_USA\_Digi\_R2