



**TECHNOPATH**  
CLINICAL DIAGNOSTICS

**Multichem<sup>®</sup> ID**  
INFECTIOUS DISEASE QC

POWERED BY



## Multichem<sup>®</sup> ID-B

Providing Instrument Specific Third-Party  
Infectious Disease Quality Controls



### Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control tubes
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

### Specifications

- ◆ 24 month closed tube stability at 2°C to 8°C.
- ◆ 90 day open tube stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-B	Positive Control	4 x 5mL	SR101B
Multichem ID-B Mini Kit	Positive Control Mini Kit	1 x 5mL	SR101MB

For further information see Instructions For Use.

INFORMATION FOR USA ONLY:  
This product is not intended for use  
in donor screening testing.

*THE* QUALITY CONTROL COMPANY

# Multichem® ID-B

## Analyte List

### Analytes

Anti-HIV 1 IgG  
 Anti-HCV IgG  
 Anti-HBc IgG  
 HBsAg  
 Anti-HTLV IgG  
 Anti-Treponema Pallidum IgG

## Analyser Systems

MANUFACTURER	SYSTEM	ASSAY NAME
Abbott	ALINITY i	Anti-HBc II
		Anti-HCV
		HBsAg Qualitative II
		HIV Ag/Ab Combo
		rHTLV I/II
		Syphilis TP
Abbott	ALINITY s	Anti-HBc
		Anti-HCV
		HBsAg
		HIV Ag/Ab Combo
		HTLV I/II
		Syphilis
Abbott	ARCHITECT	Anti-HBc II CMIA
		Anti-HCV CMIA
		HBsAg Qualitative II CMIA
		HIV Ag/Ab Combo CMIA
		rHTLV-I/II CMIA
		Syphilis TP CMIA
DiaSorin	LIAISON	Anti-HBc CLIA
		Treponema Screen CLIA
DiaSorin	LIAISON XL MUREX	HBsAg Quant CLIA
		HCV Ab CLIA
		HIV Ab/Ag CLIA
		recHTLV-I/II CLIA



**TECHNOPATH**  
 CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS  
 Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, V94 FF1P, Ireland.  
 Website: [www.technopathcd.com](http://www.technopathcd.com)