



TECHNOPATH
CLINICAL DIAGNOSTICS

COMING SOON

Multichem[®] ID-COVID19 M

Third-Party Quality Control for IgM Antibodies
to SARS-CoV-2

Features

- Ready to use liquid stable controls
- Technical and scientific support

Specifications

- 6 month* shelf life at 2°C to 8°C.
- 10 day* open vial stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-COVID19 M	Positive Control	4 x4mL	CVM200P
Multichem ID-COVID19 M	Positive Control	4 x 2mL	CVM100P

Analytes

SARS-CoV-2 IgM

For further information see Instructions For Use.

*Initial stability at launch. Will be extended once further data is available.

THE QUALITY CONTROL COMPANY

SEROLOGY
CONTROLS

Multichem[®] ID
INFECTIOUS DISEASE QC



TECHNOPATH CLINICAL DIAGNOSTICS
www.technopathcd.com
multichemcovid@technopathcd.com

MM_IDCOVID19_M_Digi_R2