

TECHNOPATH CLINICAL DIAGNOSTICS

Multichem[®] U

Providing Third-Party Test Consolidation for Urinary Chemistry QC in a Liquid Stable Format

TECHNOPATH



1



10ml TECHNOPATH, Ballina, Co. Tipperary,

Features

- Human based proprietary formulation.
- Mimics performance of patient samples.
- Targeted at key clinical decision points.
- Liquid for ease of use.

Specifications

- ◆ 24 Month closed vial stability at 2°c to 8°c.
- ♦ 30 Day open vial stability at 2°c to 8°c.

Product	Description	Configuration	Part Code
Multichem U	Unassayed Single Level (Level 1)	15 × 10mL	UC201X
	Unassayed Single Level (Level 2)	15 × 10mL	UC202X
	Assayed Single Level (Level 1)	15 × 10mL	UC201A.10
	Assayed Single Level (Level 2)	15 × 10mL	UC202A.10

THE QUALITY CONTROL COMPANY

Multichem[®] U Urine Chemistry Analyte List

General Chemistry	Pituitary/Adrenal
Amylase	Cortisol*
Calcium	Reproductive/Fertility
Chloride	
Creatinine	hCG*
Blucose	
lagnesium	
licroalbumin	
Osmolality*	
Phosphorus	
Potassium	
Sodium	
Specific Gravity*	
Jrea Nitrogen	
Jric Acid	
Jrinary Protein	



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in realtime. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathcd.com



CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS Technopath Life Sciences Park, Fort Henry, Ballina, Co.Tipperary, V94 FF1P, Ireland. For Multichem QC Support: gcsupport@technopathcd.com Website: www.technopathcd.com