



TECHNOPATH
CLINICAL DIAGNOSTICS
USA

Multichem[®] WBT

Providing Third-Party Test Consolidation for
Whole Blood Immunosuppressant QC in a
Liquid Stable Format



Features

- ◆ Human based proprietary formulation.
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 30 month closed vial stability at -20°C to -80°C.
- ◆ 10 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem WBT	Unassayed Tri-Level Multichem WBT Control	3 x 4 x 2mL	WB000X
Multichem WBT	Assayed Tri-Level Multichem WBT Control	3 x 4 x 2mL	WB000A

THE QUALITY CONTROL COMPANY

Multichem[®] WBT Immunosuppressant Analyte List

Chemistry

Folate*

Glucose*

Immunosuppressants

Cyclosporine

Sirolimus

Tacrolimus

*Please refer to lot specific package inserts for stability and performance claims.



Designed to complement and support TECHNOPATH's
Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyze their QC results in real-time. IAMQC Software products allow users to automate, centralize and standardize QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathusa.com



TECHNOPATH
CLINICAL DIAGNOSTICS
USA



TECHNOPATH CLINICAL DIAGNOSTICS USA
99 Lafayette Drive, Suite 179 Syosset, NY 11791
For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathusa.com