



TECHNOPATH

CLINICAL DIAGNOSTICS

Multichem[®] IA

Providing Third-Party Test Consolidation for Immunoassay QC in a Liquid Stable Format



Features

- ◆ Human based proprietary formulation.
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 36 month closed vial stability at -20°C to -80°C .
- ◆ 10 day open vial stability at 2°C to 8°C .

Product	Description	Configuration	Part Code
Multichem IA	Unassayed Tri-Level	3 × 4 × 5mL	IA300X
	Assayed Tri-Level	3 × 4 × 5mL	IA300A

Multichem[®] IA

Immunoassay Analyte List

Cancer Markers

Alpha Fetoprotein
Carcinogenic Embryonic Antigen
Prostate Specific Antigen, Free
Prostate Specific Antigen, Total

Cardiac

BNP
CK-MB
Myoglobin
NT-pro BNP
Troponin I
Troponin T
Ultrasensitive CRP*

Allergy

IgE

Thyroid

Anti-Thyroperoxidase
Anti-Thyroglobulin
Calcitonin
Thyroglobulin
Thyroid Stimulating Hormone
Thyroxine Binding Globulin*
Thyroxine, Free (FreeT4)
Thyroxine, (Total T4)
Triiodothyronine, Free (FreeT3)
Triiodothyronine, (Total T3)
T Uptake

Anaemia

Erythropoietin (EPO)
Ferritin
Folate
Vitamin B12

Pituitary/Adrenal

Adrenocorticotrophic hormone (ACTH)
Aldosterone*
Androstenedione*
Cortisol
Human Growth Hormone

Drugs

Acetaminophen
Amikacin
Caffeine*
Carbamazepine
Carbamazepine, Free*
Chloramphenicol*
Cyclosporine*
Digoxin
Disopyramide*
Ethosuximide*
Gentamicin
Ibuprofen*
Lidocaine*
Lithium
N-Acetyl procainamide*
Phenobarbital
Phenytoin
Phenytoin, Free*
Primidone*
Procainamide*
Quinidine*
Salicylate
Theophylline
Tobramycin
Valproic Acid
Valproic Acid, Free*
Vancomycin

Renal

Angiotensin*
Renin*

Bone Metabolism

Ostase*
Parathyroid hormone (PTH)
Procollagen NP Type 1*

Reproductive/Fertility

DHEA Sulfate
Estriol, Free
Estriol, Total*
Estrogen, Total*
Estradiol
Follicle Stimulating Hormone
Human Chorionic Gonadotropin
17-Hydroxyprogesterone*
Leutinizing Hormone
Progesterone
Prolactin
Sex Hormone Binding Globulin (SHBG)
Testosterone
Testosterone, Free*

Diabetes

C-Peptide
Insulin
Insulin-like Growth Factor (IGF-1)*

Esoterics

25 (OH) Vitamin D
Homocysteine

*Please refer to lot specific package inserts for stability and performance claims.



TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS

IRELAND
Technopath Life Sciences Park,
Fort Henry, Ballina, Co.Tipperary, V94 FF1P

USA
99 Lafayette Drive,
Suite 179 Syosset, NY 11791

For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathcd.com