



# TECHNOPATH

CLINICAL DIAGNOSTICS

## Multichem<sup>®</sup> P

Providing Third-Party Supplementary Immunoprotein QC in a Liquid Stable Format



### Features

- ◆ Human based proprietary formulation.
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

### Specifications

- ◆ 36 month closed vial stability at -20°C to -80°C.
- ◆ 14 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem P	Unassayed Single Level	12 × 3mL	SP40PX
	Assayed Single Level	12 × 3mL	SP41PA

# Multichem<sup>®</sup> P

## Immunoprotein Analyte List

### Immunoproteins

Alpha-1 Acidglycoprotein  
Alpha-1 Antitrypsin  
Alpha-2-Macroglobulin\*  
Antistreptolysin O (ASO)\*  
ADNase B (Anti-Streptococcal DNase B)\*  
Antithrombin III\*  
Apolipoprotein A1 (APO A1)  
Apolipoprotein B (APO B)  
Beta-2 Microglobulin  
C1 Inhibitor\*  
CH50 (Total hemolytic Complement)\*  
Cystatin C\*

### Chemistry Analytes

Albumin\*  
Angiotensin Converting Enzyme\*  
Total Protein\*  
Complement C3  
Complement C4  
Ceruloplasmin  
C-Reactive Protein  
Ferritin\*  
Haptoglobin  
Hemopexin\*  
Immunoglobulin A  
Immunoglobulin G  
Immunoglobulin M

### IgE\*

IgG1, Subclass\*  
IgG2, Subclass\*  
IgG3, Subclass\*  
IgG4, Subclass\*  
Kappa Light Chain\*  
Lamda Light Chain\*  
Lipoprotein (a)\*  
Prealbumin  
Properdin Factor B\*  
Retinol Binding Protein\*  
Rheumatoid Factor  
Transferrin  
sTIR (Soluble Transferrin Receptor)\*

\*Please refer to lot specific package inserts for stability and performance claims.



Designed to complement and support TECHNOPATH's  
Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit [www.technopathcd.com](http://www.technopathcd.com)



**TECHNOPATH**  
CLINICAL DIAGNOSTICS

### TECHNOPATH CLINICAL DIAGNOSTICS

#### IRELAND

Technopath Life Sciences Park,  
Fort Henry, Ballina, Co. Tipperary, V94 FF1P

#### USA

99 Lafayette Drive,  
Suite 179 Syosset, NY 11791

For Multichem QC Support: [qcsupport@technopathcd.com](mailto:qcsupport@technopathcd.com)  
Website: [www.technopathcd.com](http://www.technopathcd.com)