



TECHNOPATH

CLINICAL DIAGNOSTICS

Multichem[®] WBT

Providing Third-Party Immunosuppressant
QC in a Liquid Stable Format



Features

- ◆ Human based proprietary formulation.
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

Specifications

- ◆ 30 month closed vial stability at -20°C to -80°C.
- ◆ 10 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem WBT	Unassayed Tri-Level	3 x 4 x 2mL	WB000X
	Assayed Tri-Level	3 x 4 x 2mL	WB000A

Multichem[®] WBT Immunosuppressant Analyte List

Chemistry

Folate*
Glucose*

Immunosuppressant

Cyclosporine
Sirolimus
Tacrolimus

*Please refer to lot specific package inserts for stability and performance claims.



Designed to complement and support TECHNOPATH's
Multichem Quality Control (QC) product range.

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit www.technopathcd.com



TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS

IRELAND

Technopath Life Sciences Park,
Fort Henry, Ballina, Co. Tipperary, V94 FF1P

USA

99 Lafayette Drive,
Suite 179 Syosset, NY 11791

For Multichem QC Support: qcsupport@technopathcd.com
Website: www.technopathcd.com